




hyamira®
THE TRIGGER HASYSTEM





FILLER

BIOREVITALIZERS

**TOPICAL
ENHANCEMENT**

**HEALTH AND
BEAUTY
SUPPLEMENT**



hyamira®

Skin ageing is an evolutionary process that can be traced back to a multitude of micro-changes which occur slowly, first at cellular level and then gradually, to the extracellular matrix, involving the integrity of the tissue.

This leads to obvious wrinkles in different areas of the face, a loss of volume in the fatty tissues of the face and often radical changes to expression.

All this can be countered gradually, naturally and with the certainty of persistent action, thanks to

HYAMIRA® - THE TRIGGER HA SYSTEM AN INTEGRATED MULTIFUNCTIONAL SYSTEM





FILLER

CROSS-LINKED
Hyaluronic Acid

NATURAL RESULTS

SAFE

PREDICTABLE

A viscoelastic filler, formulated to guarantee the maximum in terms of natural-looking results.

It integrates perfectly in the skin tissues, thanks to its rheological characteristics, keeping facial volume at the correct level, both statically and dynamically.

Hyamira® is reabsorbed into the tissues with distribution kinetics depending on its infiltration depth⁽¹⁾: this guarantees **natural and predictable results**.



Ingredient	HA Reticulated	HA Reticulated	HA Reticulated
Concentration	15 mg/ml	20 mg/ml	25 mg/ml
HA Molecular weight	2400 - 4000 kDa	2400 - 4000 kDa	2400 - 4000 kDa
Cross linking agent	BDDE	BDDE	BDDE
G'	130 [Pa]	267 [Pa]	440 [Pa]
Free Residual BDDE	< 50 ppb	< 50 ppb	< 50 ppb
Directions	Treatment of periocular wrinkles	Treatment of wrinkles around the lips	Treatment of deep facial wrinkles
Injection area	Superficial-medium dermis	Medium-deep dermis	Deep hypodermis or over the periosteum
Duration of results	4-6 months	8-12 months	12-14 months
Needles	2x 30G 1/2"	2x 27G 1/2"	2x 27G 1/2"

CE 0373



FULL-FACE TREATMENT

- 3 SYRINGES OF HYAMIRA FORTE
- 4 SYRINGES OF HYAMIRA BASIC
- 1 SYRINGE OF HYAMIRA SOFT

hyamira®

FILLER

CROSS-LINKED
Hyaluronic Acid

SOFT

Treatment of periocular wrinkles

G' [Pa] 130



BASIC

Treatment of wrinkles around the lips

G' [Pa] 267



FORTE

Treatment of deep facial wrinkles

G' [Pa] 440



hyamira® BODY

FILLER

BODY RESHAPING AND CONTOURING

INCREASE TONE AND HYDRATION

NATURAL EFFECT

CROSS-LINKED
Hyaluronic Acid
10mL

Hyamira® Body is the latest generation **100% Italian Hyaluronic acid body filler**.

It allows the volumes of the body to be restored by reducing skin laxity and conveying **immediate tone** and **rapid hydration**. It reshapes large volumes harmoniously, without any surgery.



Ingredient	Monophasic cross-linked HA
Concentration	20 mg/ml
HA molecular weight	1.500 kDa
Crosslinking agent	BDDE
G' (elastic modulus)	20 Pa
BDDE residual	< 2 ppm
Injection area	Hypodermis or over the periosteum


CE 0373

hyamira® BODY

Thanks to the special formulation, it guarantees **optimal hydration** and **long-lasting firming** of depressed areas of the body. Unlike more invasive body treatments, this body filler requires a very short time for normal activity to be resumed.

APPLICATION AREAS

- GLUTES
- CALVES
- ARMS
- FOREARMS
- BODY DEPRESSIONS



hyamira®

BIOREVITALIZERS

HYALURONIC ACID

NATURAL AND PREDICTABLE RESULTS




Intradermal injections of natural hyaluronic acid for **gradual skin rejuvenation via stimulation of fibroblasts** and synthesis of new **collagen and elastin**.

Recommended for correction of soft tissue contours such as **acne and other scars**.


It combats the gradual reduction in hyaluronic acid, restoring the physiological architecture that supports skin tissue.

It improves skin hydration, compactness and elasticity in a perfectly biocompatible way.

Available in two different concentrations and in two different quantities, it guarantees focused treatment that meets all specific personal beauty requirements.



CE 0373



Results on the hydration and elasticity of the skin in 18 patients who underwent biorevitalisation treatment with Hyamira: 1 injection every 2 weeks, carried out 4 times and followed by one infiltration per month. Controls, 90 days from the last infiltration, widely confirmed the positive effects in terms of hydrated, firmer skin.

mcS: hydration parameters

Mpa: stress from pressure exerted in the point.⁽²⁾



hyamira®

16 mg/1ml and 32 mg/2ml

To prevent and treat fine surface lines, dryness and dullness on the skin. Face, neck and décolleté.

hyamira®

20 mg/1 ml and 40 mg/2ml

To recreate a lasting hydration effect on the skin, improving firmness and elasticity, and reducing fine lines. For treatment of acne scars.

PARABEN FREE

READY TO USE


NANOSOME-SONICATE-TECHNOLOGY (NST)

WHAT IS NST TECHNOLOGY?

Our laboratory use a **very high quality of Hyaluronic Acid**, parabens & chemicals free that is fractionated in low molecular weight using an ultrasound sonicator.

The low molecular weight **sonicated Hyaluronic Acid** fraction is able to penetrate deep down and, thanks to its hydrophilic nature, attracts water, boosting the skin's hydration and tone and reducing fine wrinkles.

Meanwhile, the higher molecular weight Hyaluronic Acid forms a film on the skin's outermost layer, combatting water loss and hence improving the skin's texture.



hyamira®

SKIN ACTIVATOR



Hyamira® Skin Activator is a topical serum based on **Sonicated Hyaluronic Acid plus Growth Factors** which is able to reduce wrinkles increasing the skin tone. The treatment is highly recommended to recover the dryness coming to an invasive treatment or a fractional ablation, chemical peeling, etc.

Skin Activator is encapsulated in Nanosomes, microspheres that:

- protect them from endogenous proteasis
- improve their penetration into the skin, making them real "carriers"

hyamira®

SKIN REPAIR


Hyamira® Skin Repair is a topical serum based on **Sonicated Hyaluronic Acid, Perfluorodecalin, Vitamin B5 and Heptapeptide 10** which is able to reduce the redness and the inflammation of the skins.

The treatment is highly recommended to recover the dryness coming to an invasive treatment or a fractional ablation, chemical peeling, etc.

The study conducted through the in vitro evaluation of the anti-inflammatory activity of Skin Repair on human fibroblasts, keratinocyte cells which showed a decrease of ~94% of inflammation compared to the untreated cells used as negative control.

EFFICACY STUDY by the UNIVERSITY of PAVIA





PREVENTS AND COMBATS SKIN AGEING

Synergistic with and **complementary** to any type of beauty treatment.

A complete formula. Containing anti-oxidants, it is ideal for preventing and correcting the signs of skin ageing.

Elements in the formula constantly and effectively sustain the **repair** and **renewal** processes of the connective tissue.

30

DAYS

Increases
skin hydration

60

DAYS

Increases
skin hydration,
tone and elasticity.
Reduce wrinkles.

A COMPLETE FORMULA

HYDROLYZED COLLAGEN

5000 mg

Stimulates fibroblasts to produce new collagen and elastin.

Leaves skin toned, elastic and resistant, preventing the appearance of wrinkles.

HYALURONIC ACID

100 mg

Ensures physiological hydration for the skin and boosts injected treatment results.

RESVERATROL

5 mg

Has an important anti-oxidant and anti-ageing action.

ZINC

3 mg

An essential oligoelement for skin homeostasis of the cutaneous adnexa.

It helps to protect cell membranes from the harmful effect of various stimuli, working to maintain the correct levels of hydration and promoting preservation of the pool of endogenous cell anti-oxidants.

VITAMIN B6

1650 µg

Indispensable for metabolising proteins, fats and carbohydrates, which are crucial for nourished, healthy, glowing skin.

COPPER

300 µg

Protects the cells from harmful oxidation and contributes to the wellbeing of the skin and cutaneous adnexa.

BIOTIN

50 µg

An essential vitamin for epithelial cell metabolism.

hyamira®
COLLAGEN



TAKE ONE 50 ml bottle
per day



ONLY
16 KCAL
PER
BOTTLE

10 x 50 ml bottles